

# **3D Cell Culture for Enhancement of Cell Function** Chang, R., Sun, W., Drexel University, Liu, Q., Lau, W., 3D Biotek, LLC

### Introduction

- In the complex interplay between cells and ECM, cells cultured and arranged in conventional 2D monolayer on tissue culture polystyrene surfaces fail to reflect in vivo scenario with its biological subtleties.
- In 3D culture, normal fibroblast cells can move and divide more quickly, assuming the characteristic in vivo asymmetric shape
- 3D cultures also show promise as disease models. For example, only malignant breast cells cultured in 3D carry large numbers of the specific receptors while in 2D culture, both normal and malignant breast cells have similar. high levels of the receptors
- The appropriate environment in which to culture cells in 3D presents challenges including the selection of 1) biocompatible matrix material. 2) a suitable 3D structure fabrication method, and 3) fabrication process parameters optimally characterized for enhanced biological performance

# Objective

Design and model 3D in vitro cell culture models to uncover patterns of biological phenomena that more closely mirror what happens in the in vivo microenvironment enhanced cell culture and downstream tissue engineering applications.

# Process Development





Multiple design constraints for scaffold include biophysical, anatomical, biological, transport, and manufacturability requirements

# Materials and Methods

3D InsertTM is a series of disposable 3D cell culture well inserts for use with multi-well cell culture plates. Using solid freeform fabrication technology, the scaffolds are composed of struts/filaments joined together to form a welldefined porous structure, including porosity, pore size, and surface area. Constituent materials include non-degradable or degradable polymers such as polystyrene (PS) and polycaprolactone (PCL).



#### Polystyrene

- Non-biodegradable aromatic polymer. Transparent property allows convenient
- visualization under inverted microscopy · Thermoplastic with melting temperature at about 190°C - 260°C.

#### Polycaprolactone

- biopolymers in its homo-polymeric form
- at about 58°C 60°C.

#### 3D InsertTM Design Specifications

Fabricated cylindrical scaffolds measured at 20mm in diameter and 4lavered height, strut widths of 500um, with pore sizes of 100um (1050). 300µm (3050), and 500µm (5050) and plasma surface treatment (ST). Each layer was filled with the designed scaffold pattern of a 0/90 degree orientation to generate porous structure.



#### **Biological Characterization Study Design**

- 7F2 mouse osteoblast cells (200.000 cells / scaffold seeded directly onto top surface of 3D InsertTM with Pasteur pipet
- · Control Group: 200,000 cells seeded onto 2D 12-multiwell plate
- Alamar Blue Assay (for cell proliferation) and Alkaline Phosphatase Assay (for osteogenesis) performed on 3D InsertTM at regular study time points



- A semi-crystalline aliphatic polymer.
- Slower degradation rate than most
- Thermoplastic with melting temperature

# **Results of Biological Characterization**

# 3D PCL and PS versus 2D

response of cells

· Cell function in both PCL and PS 3D InsertTM significantly better than 2D monolaver PCL material selection and smaller pore size (300µm)

confers improved osteogenic

PCL 50 PS 3050

phatase Activity as Function of Time



# 3D PS versus 2D

- · Cell function in both surface treated and non-treated PS 3D InsertTM significantly better than 2D monolayer
- Plasma surface treatment and smaller pore size (100µm) shows highest rates of proliferation and cell osteogenic activity

# Conclusions

- 3D cell culture systems, such as 3D InsertTM, mimic 3D in vivo scenario and confer improvements in cell expansion and differentiated function over the conventional 2D methods
- The precipitous decrease in cell proliferation and function over time in 2D scenario is attributed to over-confluence of cells on 2D surface where cells undergo contact inhibition and loss of viability. Also, 3D InsertTM have increased surface areas compared to 2D cell culture plates, allowing more cells to be cultured using the same size cell culture platform.
- Plasma surface treatment and smaller pore sizes confer improvements in cell proliferation and cell-specific function
- Results show reproducibly manufactured porous PS 3D InsertTM as preferred alternative for enhanced cell culture and PCL InsertTM as good candidate for enhanced cell culture and tissue engineering applications

Laboratory for Computer Aided Tissue Engineering, Drexel University http://www.mem.drexel.edu/cate

# Acknowledgement

**NSF ENG-DMI- 0427216**