

3D Perfusion Bioreactor™

User Manual and 3D Cell Perfusion Handbook





Note

This manual will be revised and updated without prior notice. To check if the latest version comes with your 3D Bioreactor System, please visit our website at http://www.3dbiotek.com/web/prod_3dbioreactor.aspx

Registration

Registration of your purchase is recommended in order to receive product update and promotional discount code for PS/PCL scaffolds specifically made for use in our 3D Bioreactor. Please follow the link below to register. <u>http://www.3dbiotek.com/web/register4bioreactor.aspx</u>





Contents

1.	3D I	Perfusion Bioreactor™ Overview4
2.	3D	Insert™ Technology5
3.	3D	Perfusion Bioreactor: Parts and Assembly6
3	3.1	3D Perfusion Bioreactor Parts List7
3	3.2	Perfusion Bioreactor Assembly8
4. (Cell Cu	ulture: Static Cell Seeding and Dynamic Cell Culture11
2	4.1	Direct Static Cell Seeding into Bioreactor Chamber12
5. F	Perfus	sion Velocity and Revolution per Minute (RPM) Graph and Table
6. F	Pump	Tubing Sterilization16
7.1	Fubin	g Connection Circuit Diagram16
8.	Арр	pendix19
8	8.1	Pre-Installation Check List19
8	3.2	Pump's Power Cord Through Incubator's Double Doors

For Latest Updates on the Bioreactor Manual, Please Visit our Website at http://www.3dbiotek.com/web/prod_3dbioreactor.aspx



1. 3D Perfusion Bioreactor™ Overview

Congratulations on the purchase of this 3D perfusion bioreactor!

This bioreactor is specifically designed for conducting physiology relevant cell culture in a 3D in vivo-like microenvironment under dynamic perfusion conditions.

This bioreactor has integrated our proven 3D Insert technology with a unique dynamic perfusion system to allow mammalian cells to grow in a 3D in vivo-like micro-environment. Users will obtain the following distinct advantages by using this 3D perfusion bioreactor.

1. Wide range of 3D surface areas for initial cell seeding

3D Perfusion Bioreactors are available in different configurations that provide cell growth area from 9.5 cm^2 to 3960 cm^2 , wide enough to meet various lab research needs.

2. Compatible with wide range of 3D Insert[™] scaffolds

A very unique feature about the 3D Perfusion Bioreactor is that a wide range of 3D Insert[™] scaffolds, from 24 well to 6 well configurations, can be used in the bioreactor, all with a simple change of a chamber.

3. Autoclavable

All parts and tubing (except pump) are autoclavable.

4. Forced CO₂ and O₂ exchange to maintain better cell growth environment

 CO_2 and O_2 are forced to bubble through the cell culture medium to better maintain an ideal cell culture environment.

5. Adjustable cell culture medium flow rate to better meet individual needs

Cell culture medium perfusion flow rate can be adjusted to meet individual needs. This also makes it possible to study the flow rate or shear force effects on cell growth. The pump's digital display and control make it easier to adjust the flow rate.



2. 3D Insert[™] Technology

One of the novelties of the 3D Bioreactor is the 3D Insert[™] porous polymer scaffolds that sit inside the bioreactor chambers. The 3D scaffold is engineered using 3D Biotek's Proprietary Precision Microfabrication Technology. Unlike other sponge like porous scaffolds, the porous 3D Insert[™] scaffolds have 100% open and interconnected pores, therefore, cell culture medium can easily perfuse through the porous structure to eliminate the poor nutrient/waste exchange phenomena often encountered with sponge-like scaffolds. The scaffolds are available in both biodegradable and non-biodegradable polymers. The entire scaffold is composed of polymer struts of uniform size. Struts on different layers of the scaffold are 90° relative to each other creating a 3-dimensional and 100% interconnected porous structure (Figure 1).

Other important features of the 3D scaffolds:

- Pre-sterilized and Ready to Use
- Well Defined Pore and Porous Structure
- Improved Cell Culture Efficiency
- Easy Separation of Cytokines and Secreted Growth Factors
- Produced from an Organic solvent Free Fabrication Process



Figure 1. CAD images of 3D InsertTM – PS scaffold. Different color struts represent different layering in the scaffold. This 4-layered model shows a top view (a) and (b) and side view (C) shows from bottom to top, layer 1 (Green), layer 2 (Yellow), layer 3 (Aquamarine), and layer 4 (Red). Each layer is 90° and offset in reference to each other.



3. 3D Perfusion Bioreactor: Parts and Assembly

Note: The 3D Perfusion Bioreactor is designed to use within a cell culture incubator at 37°C. Please do not use the pan humidifier of the incubator because the pump is not designed for use in a high humidity environment. The bioreactor comes with its own humidifier to reduce the evaporation of the cell culture medium.

The 3D Perfusion Bioreactor consists of 4 independent autoclavable polycarbonate chambers. The chambers are available in different sizes to accommodate 3D Insert^m of 24, 12, and 6-well configurations. An 8-channel peristaltic pump, which has an adjustable speed range from 0 to 60 RPM, to recirculate cell culture medium through the 4 cell culture chambers at the same time. The system is equipped with a build-in CO₂ humidifier system that constantly pumps CO₂/O₂ directly into the cell culture media.

Figure 2 and Table 1 give a detailed list of part number and description of all individual components of the 3D Perfusion Bioreactor.



Figure 2. Illustration of bioreactor Parts. Bioreactor's max dimensions Width (X) ~14", Depth (Y) ~13", and Height ~12"(Z).



3.1. 3D Perfusion Bioreactor Parts List

Table 1. Catalog number and description of 3D Perfusion Bioreactor Parts

ltem#	Part#	Description
1	3DB-10-01	Peristaltic Pump with 8 channels.
2	3DB-10-02	Aluminum Base Frame.
3	3DB-10-03	Aluminum stand for holding 4 bioreactor chambers (item# 4), 4 glass reservoirs (item# 7) & 4 small humidifier glass bottles (item# 10).
4	3DB-10-04	Bioreactor Polycarbonate Chambers (4)
5	3DB-10-05	Silicon Rubber O-rings (11/chamber x 4 chambers) to separate individual scaffolds within the chamber.
6	3DB-10-06	Pre-packagedsets of tubing (4) to close-loop the chambers with the media reservoir and Co ₂ humidifier glass bottle through the pump. Each packaged sethas 2 pre-assembled tubing for media and CO2 perfusion and 3 pieces of clear tubing (1 mid-size with white stoppers on one end and 2 short tubing pieces for the air filters. There are additional barbed tube fittings (5) for connecting bioreactor tubing circuitry during bioreactor equilibration.
7	3DB-10-07	Glass cell culture medium reservoir (4).
8	3DB-10-08	Cap with fittings (4) for glass cell culture medium reservoir.
9	3DB-10-09	Syringe Filter (0.20um) (3DB-10-07)(item# 12).
10	3DB-10-10	Co_2 Humidifier glass bottle (4). Humidified CO_2 is pumped into the medium reservoir (item# 7).
11	3DB-10-11	Cap with fittings (4) for Co ₂ Humidifier glass bottle (item# 10).
12	3DB-10-12	PS or PCL scaffolds (40) (not included) for selected chamber size.
13	3DB-10-A1	Accessory: Chamber stand (4) for holding the bioreactor chamber while loading 3D scaffolds.
14	3DB-10-A2	Accessory: Scaffold/O-ring pusher(1) for loading scaffolds and O-rings into the bioreactor chamber.
15	3DB-10-A3	Accessory: Regular forceps (1) for scaffolds and O-ring loading into the bioreactor chamber.



3.2. Perfusion Bioreactor Assembly

Note: All 3D Perfusion Bioreactor parts, <u>except peristaltic pump</u>, need to be autoclaved before any experiment. The autoclavable parts can be autoclaved using autoclaving bags at 121°C for 30 minutes and 30 minutes drying. The bioreactor tubing and syringe filters are autoclavable 5 cycles maximum.

For further technical questions on autoclavable syringe filters, visit Millipore's Technical Library at http://www.millipore.com/search.do?q=SLFG02550#0:0

The position of reservoir bottles need to be organized as 1, 2, 3, and 4 starting from left to right when assembling the 3D Perfusion Bioreactor (Please refer to Figure 2A). The perfusion bioreactor assembly has been optimized to achieve media equilibration and bioreactor chamber static cells seeding at the same time. These two steps are required for final 3D perfusion.

Each non-sterile pre-assembled tubing set comes with different lengths to account for the positions of the media reservoir and humidifier bottles on the aluminum metal stand. Make sure you don't mix all the different pre-assembled pump tubing before or after an experiment.



Figure 2A. Designation of Bottle Location

Autoclaving of all the parts, except the pump and 3D Insert[™] scaffolds, is needed before assembling the bioreactor system because the polycarbonate chambers (4), tubing, fittings, Co₂ humidifier bottles, and

media reservoir bottles are non-sterile. Before autoclaving, mark the autoclaving bags containing the tubing with the number associated to the position of the media reservoir bottle to avoid tubing mismatch. Assembly of the Perfusion Bioreactor should be done in a sterile hood. Refer to figure 2 for visual description of bioreactor tubing assembly:

- 1. Wipe the sterile hood surface with 70% ethanol prior to assembly and let the hood surface air dry.
- 2. Place all the sterile parts inside the hood. Wipe the aluminum stand holder, pump and pump's cassettes (Unlock the cassette channels from the pump by clicking downward on the oval mound on each cassette (distal end of the cassette) with 70% Ethanol. Insert the aluminum stand holder into the aluminum base frame (Refer to Figure 2 for assembly). Media reservoir bottles (item 7) should be placed in the front of the aluminum stand holder (item 3).
- 3. Place the Co₂ humidifier chambers (item 11) on the back holders in the aluminum stand. Fill each humidifier bottle with 25-ml of sterile water.
- 4. <u>3D Perfusion Tubing Bioreactor Description</u>. From this point on refer to figure 3 for assembly of media and CO₂/O₂ tubing. The tubing has been labeled with autoclavable color tape to match with the imprinted letter on the media reservoir (Red (R), Orange (O), Yellow (Y), Green (G)) and Co₂ humidifier metal caps (Green (G) and B (Blue)). Start tubing assembly with reservoir bottle 1 (from left to right). Spray the autoclave bags and bottles with 70% ethanol and wipe. Below is an itemized description of the tubing, followed by detailed description on tubing assembling.
 - Pre-assembled Bioreactor Chamber <u>Medium Loop</u>. This tubing is composed of small clear tubing, orange-labeled on one end (Orange arrow (Fig. 3)). This piece is attached to beige bioreactor tubing (Dotted black lines (Fig. 3)), which connects to large clear tubing (Black arrow (Fig. 3)) with red media stoppers on one end (Media In towards bioreactor chamber). As pre-assembled tubing, the <u>Media In towards bioreactor chamber (Red solid line, Fig. 3)</u> is connected to medium size clear tubing using a <u>large plastic tube fitting</u>. The medium size clear tubing has white media stoppers and red-labeled at the end. This will connect to the media reservoir cap (Red arrow (Fig. 3)).
 - Pre-assembled Reservoir <u>CO₂/O₂ In.</u> This tubing is composed of small blue-labeled end clear tubing connecting the humidifier cap (Blue arrow (Fig.3)). This piece is attached to beige bioreactor tubing (Dotted black lines (Fig. 3)), which connects to large green-labeled clear tubing connecting to the media reservoir cap (Green arrow (Fig. 3)).
 - \circ <u>CO₂/O₂</u>Filters. Short Clear Tubing to attach Syringe Filters (Yellow and Green solid lines with circle ends (Fig. 3)).
- 5. The following is the methodology to assemble and connect all bioreactor tubing.
- 6. <u>Bioreactor Cell Culture Media Loop</u>: connect the orange-labeled end of the tubing to the letter O (Orange) on the media reservoir cap. Secure the tubing using the left metal hook. Unhook cassette number 1, the oval mound always on the opposite side. Each cassette has an internal channel where the beige tubing will pass through. The beige pump tubing has two plastic stoppers. The first stopper secures the tubing on one end. Following, place the tubing through

the internal channel. Then, the second stopper will secure the beige tubing at the distal end. At this point, the beige tubing should be under tension. Place the cassette back into the pump, hooking the proximal side (opposite to the oval mound) first and then pressing down at the oval end until there is a clicking sound. Bring the remaining length of the tubing from the back to the front, towards the media bottle reservoir. First secure the tubing on the right hook on the metal stand. Second, pass the tubing through the bioreactor holder hoop. Then, connect the red taped end with the letter R (Red) on the media bottle reservoir. This completes the medium loop without the bioreactor chamber. Later, the bioreactor chamber will replace the <u>large plastic tube fitting connecting the Media into Bioreactor and Media out to Reservoir Bottle</u>.

- 7. <u>Bioreactor CO_2/O_2 Gas Loop</u>: connect the blue-labeled end of the tubing to the letter B (Blue) on the humidifier bottle. Unhook cassette number 2 (following cassette, next to media perfusion tubing), pressing the oval mound always on the opposite side. Each cassette has an internal channel where the beige tubing will pass through. The beige pump tubing has two plastic stoppers. The first stopper secures the tubing on one end. Following, place the tubing through the internal channel. Then, the second stopper will secure the beige tubing at the distal end. At this point, the beige tubing should be under tension. Place the cassette back into the pump, hooking the proximal side (opposite to the oval mound) first and then pressing down at the oval end until there is a clicking sound. Bring the remaining of the tubing from the back to the front, towards the media bottle reservoir. Pass the tubing through the right metal hook, connecting the green-labeled end with the letter G (Green) on the media reservoir. This completes the Bioreactor CO_2/O_2 Loop; humidified CO_2/O_2 will be actively pumped from the humidifier bottle to the media reservoir bottle during perfusion.
- 8. <u>CO₂/O₂ Filters</u>: the remaining short clear tubing, i.e. yellow and green-labeled, place a 0.20-um pore syringe filter on each end and connect the color-labeled end to the media reservoir and the humidifier bottles, yellow (Y) on media reservoir bottle and blue (B) on humidifier bottle.
- 9. Repeat steps 6 through 8 for the remaining bioreactor bottles from left to right or 1 to 4.
- 10. After connecting all tubing, fill the media reservoir bottles with cell culture media (Media bottle max. volume 90-ml). Table 2 gives approximate values of bioreactor chamber and media perfusion tubing pre-loading volumes. This take into account the extra volume to add, in addition to media for extended cell culture conditions.
- 11. Remove the silicon stoppers and place on sterile surface. Using a pipette draw desired amount of media and dispense into the media reservoir bottle. After completion, wipe the rim of the media orifice and the silicon stopper with 70% ethanol.
- 12. Place the silicon stopper tightly into the media reservoir orifice.

Table 2. Approximate volume required to fill each bioreactor chamber format and tubing

Chamber/Tubing	Total Volume
6-Well	45-ml
12-Well	15-ml
24-Well	7.0-ml
Tubing	5.0-ml



13. At this point, the bioreactor is ready for media equilibration. Place the bioreactor into a humidityfree incubator and set the pump to 20 rpm. While the media equilibrates, the next is to seed cells into the bioreactor chamber following with 3 hours static seeding.





4. Cell Culture: Static Cell Seeding and Dynamic Cell Culture

A successful 3D dynamic cell culture using 3D Perfusion Bioreactor requires the following steps:

- 1. Media Equilibration prior to 3D Perfusion
- 2. 3D Bioreactor Cell Seeding
- 3. Short Term Slow 3D Perfusion
- 4. Extreme aseptic conditions

Figure 4 is a suggested time line for conducting a 3D cell culture experiment using 3D Perfusion Bioreactor.





Figure 4. Time line for a successful perfusion bioreactor cell culture.

While the media equilibrates in the bioreactor, there are general steps to successfully set up a 3D dynamic culture cell system.

4.1. Direct Static Cell Seeding into Bioreactor Chamber

- Depending on the seeding density of the experiment, cells need to be expanded under regular tissue culture conditions. Tables 3 and 4 show the different surface areas of the 3D Insert as a function of scaffold format, fiber diameter, and fiber-to-fiber configurations. These tables help determine the seeding density on 3D scaffolds from already established 2D seeding densities into seeding densities.
- 2. Follow standard protocols for cell detachment and cell counting. Make sure the total cell suspension volume is scaled up for your specific cell number needs and 3D Biotek's suggested seeding volume. Refer to Table 5 for seeding volume values.
- 3. Once cells are in suspension, take the bioreactor chamber and seal the bottom with the provided short clear tubing with a stainless steel ball. Then, place the bioreactor chamber on the chamber stand (#13). On the opening of the bioreactor chamber cap place another small clear tubing, following with a 0.20-um air syringe filter. These will prevent any potential contamination in transit to the incubator and during the three-hour incubation period.
- 4. Using the forceps, pick an O-ring and place it on the base of the bioreactor chamber. Using the O-ring pusher, slide the O-ring evenly using rotating motion until it reaches the bottom of the bioreactor chamber.
- 5. Next, grab a 3D Insert[™] scaffold using the appropriate forceps along the scaffold's diameter and from the easy-handle sides (small cut outs for stable scaffold grip). Slide without releasing the scaffold until it evenly touches the first O-ring at the bottom. Release the forceps and retract them from the chamber.
- 6. Take another O-ring and repeat step 4.

- 7. Using the pipette, aliquot the correct volume (Table 5) containing the specified number of cells. Center the pipette tip and lower it carefully and gently until touching the scaffold (Do not push, polystyrene scaffolds will break). Release the cell suspension volume. Bring the pipette tip up.
- 8. The previous, O-ring, scaffold, O-ring, cell seeding is a sequence that ensures minimal disturbance to the seeded cells below for piece of mind (Figure 4 3D Seeding).
- 9. Repeat steps 5 through 7 until loading a total of 10 scaffolds and 11 O-rings for each chamber (Figure 4 3D Seeding).
- Close tight the bioreactor chamber using the bioreactor cap (make sure it has the syringe filter). Bring the bioreactor chamber in its holder to a humidified incubator at normal culture conditions for 3 hours static seeding.
- 11. After 3 hours, turn off the bioreactor and bring it out from the non-humidified incubator into the sterile hood. Then, bring all the 3D seeded and loaded bioreactor chambers to the sterile hood as well.
- 12. Before removing the large plastic fitting (Refer to Section 3.2, step 6), press the red and white plastic stoppers on each side of the large plastic fitting before unplugging each clear tubing end. This is to maintain constant pressure across the bioreactor.
- 13. Grab one bioreactor chamber; remove the bottom clear tubing with the stainless steel ball. Secure the bioreactor chamber on the metal stand hoop and tighten using the metal slider.
- 14. Unhook the tubing end with the red stoppers. Then, connect it to the bottom of the bioreactor chamber.
- 15. Next, remove the clear tubing with the syringe filter from the top of the bioreactor chamber. Remove the large tube fitting from the other end of the remaining tubing (with the white stoppers). Connect this tubing end to the top of the bioreactor chamber. At this point the bioreactor chamber is ready for 3D perfusion (See Figure 5 to see the large tube fitting position and bioreactor chamber for 3D perfusion (white and red asterisks)).
- 16. Repeat steps 12-15 for each bioreactor chamber.



Figure 4. 3D Seeding and 3D Perfusion Diagram



17. After connecting all 4 bioreactor chambers, unhook the red and white plastic stoppers (Warning: Failure to release the stoppers will cause pressure buildup during perfusion and experimental failure). Bring the bioreactor to a humidity-free incubator and set rpm to 1.9 to 2.1 (Figure 4, 3D Perfusion). Press start and allow slow perfusion for 24 hours. Figure 6 and table 8 give the values of velocities as a function of chamber format and rpm.



Media Equilibration

Bioreactor Perfusion





Scaffold Format	2D Growth Area	3D Growth Area	3D Growth Area	3D/2D Ratio	3D/2D Ratio
		(PS1520)	(PS3040)	(PS1520)	(PS3040)
24-Well	1.90 cm ²	10.20 cm ²	9.56 cm ²	5.4	5.0
12-Well	4.00 cm ²	21.08 cm ²	19.65 cm ²	5.3	4.9
6-Well	9.60 cm ²	54.02 cm ²	52.10 cm ²	5.6	5.4

Table 3. Total surface growth area on 3D PS scaffold

Table 4. Total surface growth area on 3D PCL scaffolds

Scaffold Format	2D Growth Area	3D Growth Area	3D Growth Area	3D/2D Ratio	3D/2D Ratio
		(PCL3030)	(PCL3050)	(PCL3030)	(PCL3050)
24-Well	1.90 cm ²	18.28 cm ²	13.74 cm ²	9.6	7.2
12-Well	4.00 cm ²	39.27 cm ²	27.90 cm ²	9.8	7.0
6-Well	9.60 cm ²	99.21 cm ²	75.62 cm ²	10.3	7.9

Table 5. Seeding volume on scaffold as a function of scaffold format

Scaffold Size	Seeding Volume (µl)					
	PS(1520)	PS(3040)	PCL(3030)	PCL(3050)		
24-well	60	120	130	160		
12-well	120	240	270	340		
6-well	300	600	690	880		





5. Perfusion Velocity and Revolution per Minute RPM Graph and Table

Figure 6. Range of rmp value for a range of biologically relevant flow velocities

In figure 6, the rpm values were calculated from flow velocity of interest and the scaffold cross-sectional area by using following equation (EQ. 1)

Where, K = 8.547, this is an experimentally determined constant from the linear fit between flow rate (ml/min) and rpm for the 3D bioreactor system ($R^2 = 0.999$).

Scaffold Format	Cross-Section cm ²
24-Well	1.13
12-Well	2.83
6-Well	7.07

Table 7.	Cross-section	Area of Sca	iffold as a f	function of	format



	Chamber Format						
	24-Well	24-Well 12-Well 6-V					
Velocity (cm/min)	rpm	rpm	rpm				
0	0	0	0				
0.2	1.93	4.84	12.08				
0.4	3.86	9.69	24.15				
0.6	5.8	14.53	36.23				
0.8	7.73	19.38	48.31				

Table 8. Relationshi	p between	pump's rp	m and biologicall	v relevant flow	velocities
		թաութ 3 ւթ	in and biologican	y i cicvant now	velocities

6. Pump Tubing Sterilization

The tubing for the pump refers to the one with 2 stops that are used with the cartridge of the pump head. The 2 stops are attached to the tubing by adhesive (Figure 7).



Figure 7. Pump Tubing with Stops

During successive autoclave procedures under high temperature, these stops may drop off the tubing due to failure of the adhesive. Given below are some recommended methods to sterilize this particular tubing.

- Gas sterilization with ethylene oxide
- Gamma Irradiation
- Autoclave at 1 BAR and 120°C for 30 minutes
- Flush ethanol and then rinse with cell media through the tubing using the pump

Consider that the tubing will mechanically degrade over time as it is repeatedly squeezed by the rollers. There comes a time when it is necessary to replace the tubing independent of cleaning considerations.



7. Tubing Connection Circuit Diagram

The following diagram (Figure 8) shows the tubing connection circuit for the system, which should be used in conjunction with the schematic diagram (Figure 3). To connect the tubing circuit for each bottle/chamber location, also refer to Figure 2A.



Figure 8. Tubing Connection Circuit for the system



8. Appendix

8.1 - Pre-Installation Check List

Before the bioreactor system is installed, it is recommended to check if the following items are available.

- Bioreactor Parts: all the items listed in Section 3.1 on Page 7.
- Scaffolds: scaffolds ordered from 3D Biotek, or your own scaffolds.
- Cells & Media: from your laboratory
- Incubators: either option (1) or (2) below
 - (1) Tall incubation with 2 independent compartments, each with standard internal space:
 18"W x 18"D x 20+"H (see Figure 9); One compartment for Media Equilibration in non-humidified condition and one for Static Cell Seeding in humidified condition.
 - (2) 1 standard-size incubator with single compartment (required minimum internal space: 16"W x 15"D x 14"H) for Media Equilibration in non-humidified condition, and 1 desktopsize incubator (no minimum internal space requirement) for Static Cell Seeding in humidified condition.
 - For the overall dimensions of the bioreactor assembly, refer to Figure 2 on Page 6.



Figure 9. Incubator with 2 independent compartments



8.2 - Pump's Power Cord Through Incubator's Double Doors

The following diagram (Figure 10) shows the setup of the pump's power cord through the double doors of the incubator.



Figure 10. Pump's Power Cord Through Incubator's double doors